



From Organisms to Genomics: From Alive to Live (cont'd)

Saria Otani, MRes., MSc., PhD

Applied Methods in Metagenomics - 23260
Denmark, 6/09/2022



Current approach



Microbial Genomics

16S rRNA gene profiling




DNA extraction
↓
16S rRNA gene specific PCR
↓
Sequencing

Metagenomics



DNA extraction
↓
Sequencing

Whole Genome Sequencing (WGS)



Cultivation
↓
DNA extraction
↓
Sequencing


Single Cell Sequencing




Cell Isolation
↓
DNA extraction
↓
Sequencing

How do we detect/study a microbiome

16S rRNA profiling



Previously on ...



End product:

A profile of all bacterial taxa in the studied niche.

Composition, who is there!


How do we detect/study a microbiome

Previously on ...




Metagenomics

sequencing-based analysis of the entire genomic content in an environmental sample.



Extract
DNA

Sequence the
entire DNA
content



End product:

Information on the diversity
and function of all organisms
in the sample

Composition, who is there!

Roles, what do they do!

Table 1 from <https://learn.inside.dtu.dk/d2l/le/content/126041/Home>

These are real 5 microbiomes.

Can you tell from this if the host of the microbiome is a chimpanzee?

Is table 1 a 16S rRNA pyrosequencing output?

What questions can you answer (what are you interested in)? (5-8 minutes)

What is your answer?

What is the most abundant bacterium (with assigned taxon)? Is it pathogenic?

Can you calculate the taxonomic assignment success? At each level.

Table 2 from <https://learn.inside.dtu.dk/d2l/le/content/126041/Home>

Is Table 2 metagenomics output? Why?

Can you tell from this if there is a chimpanzee host amongst the samples?



Table 3 from <https://learn.inside.dtu.dk/d2l/le/content/126041/Home>

What is it? Is it metagenomic output?

A project workflow in metagenomics

Plan ahead: what exactly are you looking for? design the experiment accordingly!

Sample.



DNA extraction.

4

Library preparation.

3



Sequencing platforms.

2

Analyses: *e.g.*, diagnostics

1

Plan it this way



Genomic DNA extraction – kit example

Starting material



Cell lysis





Left after the lysis:

DNA, RNA, protein

Precipitation

Clean up





Illumina Sequencing

This is commonly known
as NGS

Previously on ...



First Generation

Shotgun Sequencing

The diagram illustrates the process of shotgun sequencing. It starts with a DNA template (3' to 5') and a primer (5' to 3'). The primer is extended using four different dideoxynucleotides (ddATP, ddCTP, ddGTP, ddTTP) to create four separate DNA fragments. These fragments are then sequenced individually, resulting in four distinct peaks on a chromatogram. The peaks are labeled with their corresponding nucleotides: T, A, C, G, T.

- Sequencing by synthesis
- High accuracy
- Long read lengths
- Relatively small amount of data generated

e.g., ABI capillary sequencer (ABI)

Second Generation

Massively Parallel Sequencing

The diagram illustrates the process of massively parallel sequencing. It shows a DNA template being amplified into many small clusters on a solid support. Each cluster is then sequenced individually, resulting in a large number of short reads. The reads are then aligned to a reference sequence, and the consensus sequence is generated. The diagram also shows four circular images representing the different nucleotides: A, T, C, and G.

- Sequencing by synthesis
- Amplified templates are generated during sequencing, reducing the requirements for starting material
- High accuracy
- Short read lengths

e.g., MiSeq (Illumina), Ion Torrent (Thermo Fisher Scientific)

Third Generation

Single-molecule Sequencing

The diagram illustrates the process of single-molecule sequencing. It shows a single DNA molecule being sequenced directly, without the need for amplification. The molecule is attached to a solid support, and the sequence is determined by measuring the fluorescence of the nucleotides as they are incorporated. The diagram also shows the use of mask hairpin sequences to protect the DNA from damage. The final step is to generate a consensus sequence.

- Single-molecule templates
- Low accuracy
- Long read lengths

e.g., Single-Molecule Real-Time (SMRT) — Sequencing (Pacific Biosciences), MinION (Oxford Nanopore Technologies)

Sequencing technology selection

- **Short read technologies**

- Illumina (MiSeq, NextSeq, HiSeq, NovaSeq...)
- Ion Torrent

- **Long read technologies**

- Pacific Biosciences (PacBio)
- Oxford Nanopore Technologies (MinION)

DNA library preparation MinIon:

Same steps, yet different ingredients and handling!

MinIon


High molecular weight gDNA

Keep in the fridge.

Wide bore tips whenever handling the DNA.

Flick tubes, no vortex.


Increase all incubation times in lib. prep. protocol.



DNA library preparation head ups - Illumina:


Low sequencing output:
low DNA input - library prep failure.

Failed clustering:
Failed library prep.





**Library prep
machines
VolTRAX**

Previously on ...




Oxford Nanopore Technologies (MinION):

Long and ultra-long reads – Best high throughput available (100 Gbp) –
Portable, fast and real-time sequencing.




1. The workflow from alive material to DNA:

- Starting materials and sample selection
- DNA extraction and challenges



2. DNA preparation for sequencing:


- DNA library preparation
- Sequencing technology selection



DNA library preparation:


Collection of steps to prepare DNA to be read by a sequencing machine.

It is the same for all platforms.





Oxford Nanopore - MinKNOW

Generates FAST5 files.




Previously on ...



5' -GTGCATT-G-CTTGAGTTCTTT-CAGTT-A-3'
 |||||
 5' -CTGCATTCTCCGTGAATTCTTTCTGAATTCT-3'
 |||||
 5' -TCGC-TT-TCCGCGAATTCCTT-GAATT-C-3'

DNA library preparation:

End library product




fragmentation



end prep./repair - A tailing




adapter ligation




clean up steps


Illumina Sequencing



A. Clustering




B. High-throughput sequencing




How do we detect/study a microbiome

Current approach



Microbial Genomics





Single cell sequencing:

Microfluidic.

Established in Illumina.


Enables targeted sequencing of genetic content within one cells.

Traceability!




1. The workflow from alive material to DNA:

- Starting materials and sample selection
- DNA extraction and challenges



2. DNA preparation for sequencing:

- DNA library preparation
- Sequencing technology selection



Illumina vs. Oxford Nanopore

– Illumina (MiSeq, NextSeq, HiSeq, NovaSeq...):

Average 300 bp reads - short

Good accuracy

Error rate ~0.1%


– Oxford Nanopore Technologies (MinION):

Very long reads (up to 900 kb.)

Fast turnaround time (down to 2 hrs)

Portable and real-time sequencing


Large error rates (<1% - 2022)




Illumina vs. Oxford Nanopore: DNA input

Illumina: DNA concentration ≈ 1 ng


Nanopore: DNA concentration > 1000 ng




Library preparation




A *Klebsiella pneumoniae*



B *Acinetobacter baumannii*



C



Example


What are the plasmidome dynamics in environmental samples like soil and sewage water?

- a. What is the first thing to do?
- b. Which DNA should be extracted?
- c. Which sequencing technology do you recommend?

Global sewage plasmidome characterisation (Nanopore)

Plasmidome dynamics in 24 sewage samples from around the globe

Sewage




High quality DNA from sewage is problematic due to sample nature:

- Low DNA content
- Mixed microbial community with extra amount of **contaminants** (block DNA extractions)

Global sewage plasmidome characterisation (Nanopore)

Pre-sequencing treatment:




Global sewage plasmidome characterisation

Input QC: DNA quality, purity and length

Ensure good quality DNA input

Goes to sequencing


DNA modified protocol (beading, longer lysis, longer protein precipitation...)



concentration
DNA purity




DNA standard protocol




Global sewage plasmidome characterisation

Comparison example, differences in DNA input between 2 sequencing technologies:


Illumina DNA concentration ≈ 1 ng
:
Nanopore: DNA concentration > 1000 ng



Library preparation



Plasmidome characterisation



Comparing plasmidome
Nanopore throughputs

Plasmidome dynamics in
24 sewage samples from
around the globe


Which protocol is the best
at exploiting long read
seq.?



**Which one is the best
output to exploit long
read sequencing techn.
the best way**


① Start presenting to display the poll results on this slide.

Global sewage plasmidome characterisation




Global sewage plasmidome characterisation

Plasmid content
in 24
samples from
22 countries
show diluted
signal of
geographical
differences.




Global sewage plasmidome characterisation (Nanopore)




Philipp Kirstahler



Frederik Teudt



Global sewage plasmidome functional characterisation



MGEfinder



- **Input:** Assembled sequence data (contigs or scaffolds)*
- Detects MGEs based on sequence similarity to known elements
- Includes a curated database of ~4,400 MGEs
- **Read more:** <https://doi.org/10.1093/jac/dkaa390>

Webserver,

- Graphical user interface
- Annotates AMR genes, Virulence factors & Plasmid re
- User friendly but less flexible



Local installation

- Python package, hosted on PyPI, installed with pip
- Flexible with user customizable thresholds



Plasmid-host prediction from sewage plasmidomes



Plasmid-host prediction using machine learning from plasmidome

“PlasmidHostFinder: Prediction of plasmid hosts using random forest”



Frederik
Teudt



Derya
Aytan-Aktug

Center for Genomic Epidemiology

Username
Password

HomeServicesInstructionsOutput

PlasmidHostFinder 1.0

Plasmid host range prediction
View the [version history](#) of this server.

Please note that the program only works with assemblies (.fasta/.fna)!

Compressed files are also not acceptable.

Select the taxonomic level that the prediction will be reported

Select the class probability threshold for the prediction

Select mode of the program

Choose File(s)

Name	Size	Progress	Status
<div><div>Upload</div><div>Remove</div></div>			

Case 3:

Patient with gastrointestinal infection and diarrhoea: Sample: Faeces. Went to the hospital and was admitted to ICU immediately.

2: If we are to do metagenomics sequencing on the faecal sample of the patient, which technology would you use if the results are needed within one day with real-time result observation:

A: Illumina

B: 454


C: Oxford Nanopore

D: PacBio

Antimicrobial resistance and bacterial pathogen detection on-site

DNA extraction and library preparation (wet lab):

Can it be done in the field?



Antimicrobial resistance and bacterial pathogen detection on-site



FARMED:

Fast detection of pathogens and AMR genes on-site



```
#!/bin/bash

# Download software, reference genome, and FASTQ files
wget http://downloads.sourceforge.net/project/bio-bwa/bwa/bwa-0.7.12_x64-linux.tar.bz2
wget -O refGenomeFiles.zip https://ndownloader.figshare.com/files/4841889
wget -O FASTQ.zip https://ndownloader.figshare.com/articles/3114454/versions/1

# Extract software and reference genome files
tar -jxvf bwa/bwa-0.7.12_x64-linux.tar.bz2
unzip refGenomeFiles.zip
unzip FASTQ.zip

# Align FASTQ data to reference genome and save to SAM file
bwa.kit/bwa mem KJ668346.fa SRR1972917_1.fastq SRR1972917_1.fastq > SRR1972917_aligned.sam

# Convert SAM file to BAM file
bwa.kit/samtools view -bS SRR1972917_aligned.sam > SRR1972917_aligned.bam


# Sort BAM file
bwa.kit/samtools sort SRR1972917_aligned.bam SRR1972917_aligned_sorted

# Index BAM file
bwa.kit/samtools index SRR1972917_aligned_sorted.bam
```

November - May

DNA extraction and sequencing in field

Sequence analyses in field



Not only in Europe

In Tanzania (different project)



Your sequence quality is of your DNA.

The better your DNA quality, the better your sequencing.

DNA extraction is material-dependent.

Know your starting material and downstream application.


Why are you sequencing.

- Short read technologies


- Illumina (MiSeq, NextSeq, HiSeq, NovaSeq...)
- Ion Torrent

- Long read technologies



- Pacific Biosciences (PacBio)
- Oxford Nanopore Technologies (MinION)



Choose your sequencing technology that serves your analysis best (*e.g.*, gene detection, identification).



To remember ☺





Next-Generation sequencing as a diagnostic tool (Microbiology)

We have come a long way



Antimicrobial Susceptibility Test



"Now, here, you see, it takes all the running you can do just to keep in the same place. If you want to get somewhere else, you must run at least twice as fast!"

- Red Queen from *Through the Looking Glass*



Thank you ☺

Twitter: @SariaOtani

LinkedIn: Saria Otani

Email: saot@food.dtu.dk

Questions 😊